

ALPS™ SMART SEAL LINER

SAFR / SADT / SBFR



Transfemoral (**SA**) rings begin 14 cm from the center of the distal end.

Transtibial (**SB**) rings begin 9 cm from the center of the distal end.

Gel Liner with Bands – With or Without Distal Attachment

The ALPS™ Smart Seal Liner features raised bands that grip the socket wall to form a secure interface. The low modulus bands stretch against the socket wall, while the smooth inner wall conforms easily to the residual limb. The Smart Seal Liner has a new seamless knitted construction with a modern look and is available in locking and cushion suspension. The locking version features a new distal construction to control distal distraction.



Locking



Cushion



High Density (HD) Gel

- New seamless knitted construction
- Only available in HD Gel
- Available in locking and cushion suspension in 3mm or 6mm uniform thickness




Covered by one or more of the following US Patent No.
6,454,812 7,709,570



Tel: +420 37 72 23 127
www.easyliner.eu | info@easyliner.eu


Cross Section Diagram

3mm Uniform, 6mm Uniform





How to Size

Measurements taken 6cm from the distal end




Sizing Chart

Transfemoral Smart Seal Suction Gel Liner

Measured Size (cm) Description	 SAFR Cushion	 SADT Locking
	3mm or 6mm Uniform	3mm or 6mm Uniform
20-23 cm		SADT20-(3 or 6)
24-25 cm		SADT24-(3 or 6)
26-27 cm	SAFR26-(3 or 6)	SADT26-(3 or 6)
28-31 cm	SAFR28-(3 or 6)	SADT28-(3 or 6)
32-37 cm	SAFR32-(3 or 6)	SADT32-(3 or 6)
38-43 cm	SAFR38-(3 or 6)	SADT38-(3 or 6)

Transtibial Smart Seal Suction Gel Liner

Measured Size (cm) Description	 SBFR Cushion
	3mm or 6mm Uniform
20-23 cm	SBFR20-(3 or 6)
24-25 cm	SBFR24-(3 or 6)
26-27 cm	SBFR26-(3 or 6)
28-31 cm	SBFR28-(3 or 6)

ALPS SOUTH EUROPE S.R.O.
Božkovské náměstí 17/21
32600 Plzeň – Czech Republic
Tel. +420 37 72 23 127 Email: info@easyliner.eu